


Marker Band Positioning & Swaging Machines



Automatic Band Positioner & Pre-Crimper Model PX

- ✓ Extremely accurate & productive
- ✓ Suitable for High-Volume Output

Simple Band Positioner & Pre-Crimper Model SP

- ✓ Quick & easy manual operation
- ✓ Lowest cost solution

Band Swaging Machine Model SG

- ✓ Most Consistent Swaged Band Diameter: no force or pressure setting = *lowest diameter variability*
- ✓ Very simple to set up & use; No shimming required
- ✓ Automatic Product Feed
- ✓ Clean-room ready
- ✓ Tray Loader for High-Volume Output (Optional)